

NADIR® Flat Sheet Membrane

Format Size Options

MICRODYN-NADIR offers NADIR® flat sheet membrane in rolls or sample sizes for testing purposes. NADIR flat sheet membrane is available in rolls of 150 m x 1016 mm (492 ft x 40.0 inches) or in the various format sizes listed below. For more information, please refer to product spec sheets or contact MICRODYN-NADIR Technical Service.

NADIR® FLAT SHEET MEMBRANE FORMAT SIZES

Format	Format Dimensions (mm)	Format Dimensions (inches)	Available Membranes
A4	200 x 300 mm	8.3 x 11.7 inches	All NADIR membranes
M20N	Outer diameter 250 mm, inner diameter 59 mm	Outer diameter 9.8 inches, inner diameter 2.3 inches	All NADIR membranes
M30N	300.3 x 119.62 mm	11.8 x 4.7 inches	UPO10, UPO20, UHO04, UHO30, UHO50
M36N	424 x 313.05 mm, 2 holes of 73 mm diameter	16.7 x 12.3 inches, 2 holes of 2.9 inches diameter	UPO10, UPO20, UHO04, UHO30, UHO50
M37N	380.12 x 308.15 mm, 3 holes of 72.86 mm / 72.79 mm / 158.51 x 86.57 mm	15.0 x 12.1 inches, 3 holes of 2.9 inches / 2.9 inches / 6.2 x 3.4 inches	UPO10, UPO20, UHO04, UHO30, UHO50
M38N	424 x 313.05 mm, 2 holes of 98 mm diameter	16.7 x 12.3 inches, 2 holes of 3.9 inches	UPO10, UPO20, UHO04, UHO30, UHO50
M39N	809.5 x 202.50 mm	31.9 x 8.0 inches	All NADIR membranes

Contact

Europe

Germany: +49 611 962 6001
Italy: +39 0721 1796201
info@microdyn-nadir.com

Americas

USA: +1 805 964 8003
sales.mnus@microdyn-nadir.com

Asia

Singapore: +65 6457 7533
China: +86 10 8413 9860
waterchina@mann-hummel.com

FORM NO.: TB-042
REVISION DATE: 02/05/2020